

## pUC57


### V1400

质粒 -20℃保存  
菌株 4℃保存

载体: pUC57

菌株: DH5 α

抗性: Amp



M13F(-47) primer → 396 EcoRI SacI KpnI Bsp68I Mph103I XbaI EcoRV BamHI SmaI ApaI Sall PstI Eco147I PaeI HindIII 476

5' G TAA AAC GAC GGC CAG TGA ATT CGA GCT CGG TAC CTC GCG AAT GCA TCT AGA TAT CGG ATC CCG GGC CCG TCG ACT GCA GAG GCC TGC ATG CAA GCT TGG  
3' c ATT TTG CTG CCG GTC ACT TAA GCT CGA GCC ATG GAG CGC TTA CGT AGA TCT ATA GCC TAG GGC CCG GGC AGC TGA CGT CTC CGG ACG TAC GTT CGA Acc  
LacZ ← Val Val Ala Leu Ser Asn Ser Ser Pro Val Glu Arg Ile Cys Arg Ser Ile Pro Asp Arg Ala Arg Arg Ser Cys Leu Gly Ala His Leu Ser Pro

CGT AAT CAT GGT CAT AGC TGT TTC CTG 3'  
GCA TTA GTA CCA G7A TCG ACA AAG GAC 5'  
Thr Ile Met Thr Met

← M13R(-48) primer

Cloning Site